

# Your safety is our priority

## Technical Specification Sheet

### KN95 Mask

#### Key Features

- Double filter
- Adjustable nose clip
- Flat-fold style
- High efficiency level of at least 95% against non-oil aerosols
- Lightweight, ear-loop mask is easy to put on or remove



#### Material Composition

- Straps - Polyisoprene
- Nose Clip - Aluminum
- Coverweb - Polypropylene
- Filter - Meltblown fabric

Not made with natural rubber latex

#### Intended Uses

- Intended to reduce wearer's respiratory exposure to certain airborne contaminants.
- Always follow user instructions and use as directed

#### Do Not Use For

- DO NOT use for gases or vapors (i.e. anesthetic gases such as isoflurane or vapors from sterilants)
- DO NOT use in any manner not indicated in the user instructions
- DO NOT use in any surgical setting or where significant exposure to liquid, bodily or other hazardous fluids, may be expected;
- Do NOT use in a clinical setting where the infection risk level through inhalation exposure is high
- DO NOT use in the presence of a high intensity heat source or flammable gas.

#### Approvals and Standards

- Meets GBT 32610-2016 requirements for a minimum 95% filtration efficiency against solid and liquid aerosols that do not contain oil.
- EN-149:2001:FFP2
- CE Certificate No: P6301727

### Ordering Information

Description	Mask/Box	Box/Package	Total Pieces	Case Size
KN95 Mask	1000	1	1000	53cm*47cm*53cm

## Time Use Limitation and Storage

- May be used until damaged, breathing becomes difficult, or contaminated with blood or body fluids.
- Store in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and chemicals
- Do NOT use for more than 20 hours
- Store in temperatures between -4°F (-20°C) and +86°F (+30°C), that do not exceed 80% RH

## Warning!

This helps protect against certain airborne contaminants but does not eliminate exposure to or the risk of contracting any disease or infection. **Misuse may result in sickness or death.** Before use, the wearer must read and understand the User Instructions.